



Disposable Protective Medical Coverall

Part No: 0100-GS GN10016



Material : PP + PE (Non-woven fabric - polypropylene + polyethylene breathable film)
Fabric weight : 62 gsm
Certifications : FDA and CE (14126:2003 / AC:2004)

Sizing Information					
Size	Height (cm)	Bust (cm)	Sleeve length	Cuff width	Ankle width
XS	165	120	84	18 cm	24 cm
S	169	125	86		
M	173	130	90		
L	178	135	93		
XL	181	140	96		
XXL	188	145	99		

Packaging Information

1pc/bag, 45 pieces/carton
Carton size : 70cm*46cm*30cm
Weight of carton : 15kg

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season Building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2