



## Public Access Bleeding Control Semi-Recessed Metal Wall Case



North American Rescue's Public Access Bleeding Control Semi-recessed Metal Wall Case provides low profile, high visibility access to critical medical supplies for use by initial first responders and bystanders in public areas. It is an ideal solution to safely store & protect your Public Access Bleeding Control Kits.

The white powder painted corrosive resistant galvanized steel cabinet is ADA compliant and features a stainless steel continuous hinge & chrome plated handle designed to open with one quick pull, without the need to twist and pull. The semi-recessed design mounts permanently inside the wall with a low profile 1 in. trim. Features include highly visible "Bleeding Control Kit" lettering for fast identification, convenient Instructions for Use posted on the interior back panel and clear plexiglass window for viewing station contents. Alarm is available upon special request. Contact Customer Service for details.

### Key features at glance

- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge & chrome plated handle designed to safely store the Public Access Bleeding Control Pack in public areas
- Provides both physical and environmental protection for medical equipment
- Semi-recessed design mounts in the wall with a low profile 1 in. trim
- ADA compliant and designed for quick one pull to open - No twist and pull

---

### Golden Season Pte Ltd



- Features highly visible “Bleeding Control Kit” lettering for fast identification and convenient Instructions for Use on the interior back panel
- Alarm available upon special request. Contact Customer Service for details

| Specifications        |   |
|-----------------------|---|
| Dimensions (Outside): | H 17.5 in. x W 17.5 in. x D 1 in.         |
| Dimensions (Inside):  | H 14.125 in. x W 14.125 in. x D 6.875 in. |
| Weight:               | 10 lb                                     |

\*Specifications are subject to change without prior notice.