



Contactless Thermometer



The contactless thermometer can be used to measure the body temperature, which is suitable for people of all ages. It can be used to measure the temperature of milk and coffee as well.

With the unique technology, the contactless thermometer can give the result with non-contact, stable measurement and without the interference from ambient temperature. The product will self-test every time it starts up to make sure the normal operation and accuracy.

Product Enhancements:

- ✓ Infrared technology specially designed to take the forehead temperature of person.
- ✓ Specially designed to take the milk temperature
- ✓ Allows for continuous temperature taking
- ✓ Capable of storing 20 sets of measurement data
- ✓ No contact hence no need to disinfect after use (no probe to change)
- ✓ Bi-color Back-lighted LCD Screen
- ✓ Providing two temperature modes °F/°C
- ✓ Fever alarm function for forehead temperature test

Golden Season Pte Ltd



Specifications

LCD display resolution (definition):	0.1
Body temperature:	32.0°C _ 42.2°C
Milk temperature:	0.0°C _ 100.0°C
Measurement error:	± 0.3°C
Repeatability:	± 0.2°C
Measurement time:	0.5s
Measurement distance:	10mm _ 60mm
Auto-shutdown time:	10s

*Specifications are subject to change without prior notice.



Golden Season Pte Ltd