



## CPR Board



The CPR Board, AKA Cardiac Board / Cardiac Arrest Board, aids in the administration of CPR by providing a flat rigid surface for CPR administration.

The product is designed to position the patient for effective application of CPR.

The CPR Board is manufactured from high density polypropylene, and is able to float on water. The item is suitable for ambulance equipping, safety stations, swimming pools etc.

Specifications	
Length:	60 cm
Width:	72 cm
Height:	7 cm
Weight:	1.35 kg
Material Type:	High density polypropylene

\*Specifications are subject to change without prior notice.

### Golden Season Pte Ltd