




## Floating Spine Board



The Floating Spine Board is buoyant to 95kgs making it suitable for water rescue operations. It is X-Ray translucent, lightweight, waterproof and resilient to wear and tear. The Floating Spine Board is compatible with our head immobilizers to secure the patient head down for casualty evacuation. In doing so, further trauma caused by head / back movement is prevented.

Specifications	
Length:	184 cm
Width:	46 cm
Height:	6.5 cm
Weight:	8.54 kg
Material Type:	High density polypropylene
Accessories:	3 patient restraints included

\*Specifications are subject to change without prior notice.

---

### Golden Season Pte Ltd