



GS5 Ambulance Stretcher



The GS5 Ambulance Stretcher is designed for city conditions where patients are commonly found on high rise apartments.

The main feature of the product is it's ability to fold down into a wheelchair, allowing for access in small elevators. The conversion ferom Ambulance stretcher to wheelchair can be done even when the patient is on the stretcher.

The GS5 is also equipped with a gas spring to allow varied resting angles for the patient. It can be operated even when the patient is on the stretcher, preventing unnecessary trauma associated with movement with injuries.

The product comes complete with IV drip stand, locking bracket, resting cushions, lockable wheels.

The legs of the stretcher fold down to allow entry into ambulances and locks with the bracket provided.

Golden Season Pte Ltd



Key features at glance:

- Able to convert to wheelchair
- Rugged and Durable
- Gas spring operated backrest
- Collapsible side rails
- Easy to operate

GS5 Ambulance Stretcher – Specifications



Specifications:

Length:	198 cm
Width:	55 cm
Height:	86 cm
Low position height:	25 cm
Resting angle:	0° to 85°
Maximum Recommended Load:	159 kg

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd