



## Kendrick Extraction Device (K.E.D.)



The Kendrick Extrication Device (KED) is used in vehicle extrication to remove victims from collided motor vehicles. It is carried on rescue vehicles, first aid responders as well as ambulance services. The KED is normally used together with our cervical collars. The aim of the KED and cervical collar is to secure the head, neck and torso of patient preventing further trauma as well as aggravating injuries during extrication.

| Specifications      |  |
|---------------------|--|
| Dimension (Deploy): | 82cm x 83cm  |
| Dimension (Folded): | 86.5cm x 39cm x 11cm   |
| Gross Weight        | 3.37 kg  |
| Material Type:      | Vinyl coated nylon   |
| Accessories:        | Wrap around vest, Adjustable neck roll, 2 head straps, carrying case |

\*Specifications are subject to change without prior notice.