



## Spinal Board



The Spinal Board is used by EMS professionals, paramedics and rescue teams to facilitate extrication of patients sustaining spine injuries or other injuries.

### Reduces trauma to patient

The Spinal Board has flat surface and does not aggravate cervical spinal injuries or compromises CPR efforts.

### Ergonomic Handles

Golden Season's Spinal Board comes with in built ergonomic handles to facilitate with lifting and patients evacuation.

### 100% X-ray, CT, and MR compatibility

Golden Season's Spinal Board is 100% X-ray, CT and MR compatibility. It is also designed to with stand abuse and rugged handling in the field.

### Specifications

Dimension:	184cm x 45cm x 6cm
Weight:	7.65kg
Material Type:	High density polypropylene
Accessories:	3 Patient restraints

\*Specifications are subject to change without prior notice.

### Golden Season Pte Ltd