



T3 Lightweight Assault Litter



Our Newest Evolution of Evacuation Platforms... Weighing less than 12 lbs!

T3™ Lightweight Assault Litter is designed to support Special Operation Forces where every pound matters. State of the art materials and technologies have allowed us to reduce the weight without compromising the overall strength and performance of the evacuation platform. The overall length of 81 in. enables floor loading on a UH-60 helicopter while allowing for maximum maneuverability during casualty evacuation in urban terrain. Weighing in at less than twelve pounds, the T3™ is the lightest rigid litter in the world.

Dimensions	
Folded:	H 21 in. x W 8.75 in. x D 7.25 in.
Open:	L 81 in. x W 22.5 in. x D 6 in.
Weight:	11 lb 13 oz (5.34 kg)
Payload:	1200 lb (544.3 kg)
Ground Clearance:	1.5 in.

Special Features:

- ✓ Weighs less than 12 lbs.
- ✓ Maximum maneuverability during TACEVAC in urban terrains
- ✓ Overall length of 81 in. enables floor loading on a UH-60 helicopter
- ✓ 2 straps for securing patients
- ✓ 4 IV attachment points
- ✓ Ergonomically designed handles
- ✓ Quad-fold design for maximum cube space storage

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd