



Talon II 90C Quad Fold Stretcher



When the U.S. Army envisioned an improved litter, North American Rescue built it. Talon® II 90C, the most compact & versatile litter on the market, allows casualty transport in restricted compartments. By extending the ergonomically designed handles, the Talon® becomes a standard NATO-compatible evacuation platform. Already the U.S. Marine Corps official battlefield litter, the Talon® 90C is today's ultimate TACEVAC solution.

Dimensions	
Folded with handles collapsed:	H 20.5 in. x W 8.5 in. x D 9.5 in.
Open with handles collapsed:	L 78 in. x W 22.5 in. x D 6 in.
Open with handles extended:	L 90 in. x W 22.5 in. x D 6 in.
Weight:	16 lb 6 oz (7.4 kg)
Payload:	1200 lb (544.3 kg)
Ground Clearance:	1.5 in.

Special Features:

- ✓ Certified Safe-to-Fly U.S. Army Rotary Wing
- ✓ Retractable, ergonomically designed handles
- ✓ Automatic locking hinges
- ✓ Slip-resistant, chemical-resistant ripstop fabric tested & approved by Federal & private labs using live warfare agents
- ✓ 6 IV-attachment points
- ✓ 2 patient-securing straps
- ✓ U.S. Patents 6,842,923 & 5,598,592, Patent(s) Pending

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd